



Implant Checklist

- Detailed work order including:**
 - Teeth to be restored with implant type, size and manufacturer. Include oral surgeon's print-out of implant/s placed.
 - Material choice for abutment and crown: Titanium, zirconia, lithium disilicate. PFM.
 - Style of restoration: screw-retained, screwmentable, UCLA castable, cementable.
 - Shade selected
 - Patient goals
- Impression of maxillary and mandibular arches**
 - Accurate impression of implant position and soft tissue contours.
 - Include impression post and lab analog.
- Pre-op models**
- Bite record**
- Model or impression of approved provisionals**, when necessary
- Photographs taken with patient sitting upright or standing**, suggested views:
 1. 1:10--Full face with natural smile, head in vertical position and nose in center of image.
 2. 1:2--full natural smile, non-retracted.
 3. 1:2--retracted, slightly parted, frontal, left and right views.
 4. 1:1 showing any surface texture, incisal details, teeth parted, slightly underexposed.
 5. Shade tab held up to the incisal edge and in the same plane as the tooth, with the shade tab label clearly visible.
 6. Any other views which are important to the successful completion of the case.

Photographs can be sent in several ways:

1. Compressed and attached to an email.
2. Burned onto a CD or loaded onto a memory stick. These will be returned with the completed case.
3. Post photos to a share photo site, such as Dropbox.